

180° ROBINSON GONIOMETER 7" - stainless steel



Product Code	24902
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279249026
---------	---------------

Description

180° ROBINSON GONIOMETER 7" - 17.5 cm

Pocket size 180° Robinson goniometer for doctor and therapist.
Linear scale marks 5° increments.

Multilanguage manual: GB, FR, IT, ES.

Approximate measurement of values.