

360° GONIOMETER 14" - stainless steel



Product Code	24903
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279249033
---------	---------------

Description

360° GONIOMETER 14" - 35 cm

This goniometer has two opposing 180° scale and one 360° scale all marked in 1° increments.
Sturdy adjustable knob, controlling tension of arms.

Multilanguage manual: GB, FR, IT, ES.

Approximate measurement of values.