

# 180° EXTRA LONG ARM GONIOMETER 14" - stainless steel



Product Code	24904
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279249040
---------	---------------

## Description

180° EXTRA LONG ARMS GONIOMETER 14" - 35 cm

Extra long arms, useful for measuring all joints.  
Arm has a linear scale marked in 1° increments.

Multilanguage manual: GB, FR, IT, ES.

Approximate measurement of values.