

# WUNDER BABY 02-1 PROFESSIONAL ELECTRONIC SCALE 6-15 kg - Class III



Product Code	24990
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	2497343
CND	V0399
EAN/UPC	8023279249903
GMDN	35324

## Description

WUNDER BABY02-1 PROFESSIONAL ELECTRONIC SCALE 6-15 kg - Class III

Professional and portable electronic scale for babies, made of white ABS with digital indicator.

Functions: weight, tare, milk intake function, hold, auto zero, programmable auto-off.

Breast-milkintake function determines how much milk the baby has suckled.

Comfortable and safe surface thanks to the curved edges. Easy to read double display. MDR certified.

Supplied with power supply, RS232 interface, level, adjustable feet and safety brackets.

User manual: GB, IT, ES. Available on [gimaitaly.com](http://gimaitaly.com): FR, PT, DE, NL, SE, NO, FI.

Made in Italy.

## Technical Specifications

Capacity - sensitivity: 6-15 kg - 2/5 g

Display: LCD 20mm weight with 5 active digits / LCD 12mm milk intake

Weighing surface: 560x290 mm

Column height: 460 mm

Size - weight: 560x450x460 mm - 4,8 kg

Power supply (included): adapter 12V - 2 A

Power supply (optional): rechargeable battery or 6xAAA alkaline batteries

Cart size: 480x640x910 mm

Cart load capacity: 100 kg