

WUNDER WH80 MANUAL HEIGHT ROD 35-80 cm with analogue reading



Product Code	24992
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2703560
CND	V9099
EAN/UPC	8023279249927

Description

Manual height rod with analogue reading 35-80 cm - optional

Size 87x30xh 7 cm

Telescopic, aluminium height rod with double ABS paddle (fixed headrest and sliding footrest)