


WUNDER RB200 PROFESSIONAL DIGITAL SCALE 250 kg - Class III - standard



Product Code	25000
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	2497352
CND	V0399
EAN/UPC	8023279250008
GMDN	38211

Description

WUNDER RB200 DIGITAL SCALE - 250 kg - Class III - standard

Professional digital column scale with a large and stable painted aluminium base with non-slip painted metal platform. The scale is equipped with a dual LCD display with triple simultaneous readout: weight, height and BMI. Adjustable display for clearer and more comfortable reading. Functions: weight, tare, mother/child weight, hold, BMI, autozero, programmable auto-off. 4 large adjustable feet and 2 integrated wheels for easy transport.

Supplied with power adaptor, RS232 interface, level, 4 adjustable feet, 2 wheels and multilingual manual: GB, IT, ES, on request FR, DE.
MDR certified.

Made in Italy.

Technical Specifications

RB200 DIGITAL SCALE - 250 kg - Class III - standard

Capacity: 150/250 kg

Sensitivity: 50/100 g

Display: double LCD: LCD 20mm weight with 5 active digits; LCD 12 mm height and BMI

Scale size - weight: 275x430xh 900mm - 7 kg

Packing size: 330x1100x290 mm

Platform size: 275x295xh 55mm

Column size: Ø60 x h 800mm

Power supply (included): adapter 12V - 2 A

Power supply (optional): rechargeable battery or 6xAAA alkaline batteries