

# WUNDER DE20 ELECTRONIC CHAIR SCALE 250 kg - Class III - red - standard



Product Code	25020
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	2496387
CND	V0399
EAN/UPC	8023279250206
GMDN	35996

## Description

DE20 ELECTRONIC CHAIR SCALE - 250 kg - Class III - red - standard

High precision electronic chair scale with an innovative design.

Functions: weight, tare, hold, BMI, auto zero, programmable auto-off.

Painted metal tube structure with adjustable ABS display, 4 360° swivelling castors with 4 safety brakes, padded handle for comfortable transport.

Adjustable footrests and reclinable armrests up to 90° to facilitate patients access.

Dual LCD display with triple simultaneous reading: weight, height and BMI.

Compact packaging reduces shipping costs.

Supplied with power adaptor, RS232 interface, level, 4 castors with brakes, adjustable armrests and footrests, multilingual manual: GB, IT, ES. On request FR, DE. MDR Certified.

Made in Italy.

## Technical Specifications

Capacity: 250 kg

Sensitivity: 50/100 g

Display: Double LCD: LCD 20mm weight with 5 digits; LCD 12mm height and BMI

Size - weight: 590x1040xh 960 mm - 22 kg

Seating size: 550x400 mm

Packaging size: 620x810xh 620 mm

Power supply: 12V - 2A power adaptor (included), 6xAAA alkaline batteries or rechargeable battery (optional)