

# WUNDER RW2.0 DIGITAL PLATFORM SCALE WITH HANDRAIL 300 kg - Class III



Product Code	25026
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No
RDM (NSIS)	2497355
CND	V0399
EAN/UPC	8023279250268
GMDN	15053

## Description

WUNDER RW2.0 DIGITAL PLATFORM SCALE with handrail - 300 kg - Class III

Large and easily accessible platform with double ramp, ideal for wheelchair patients.  
Functions: weight, tare, hold, BMI, auto zero, programmable auto-off.  
Dual ABS display with triple simultaneous readout: weight, height, BMI.

Supplied with: power adaptor, RS232 interface, level, adjustable feet, double ramp, foldable handrail, carrying handle, 2 castors.  
MDR certified.

Multilingual manual: GB, IT, ES

Made in Italy.

## Technical Specifications

Capacity - Sensitivity: 300 kg - 100 g

Display: Double LCD: LCD 20mm weight with 5 digits; LCD height and BMI 12mm

Size: 900x1150xh 1100mm

Weight: 30.5 kg

Power supply: adapter 12V - 2A (included), 6xAAA batteries (not included) or rechargeable battery pack (optional)

## Standard accessories

<div>Power adaptor</div><div>RS232 interface</div><div>Level</div><div>Adjustable feet</div><div>Double ramp</div><div>Carrying handle</div><div>Foldable handrail</div><div>Manual: GB, IT, ES.</div><div>Available on gimaitaly.com: FR, PT, DE, NL, SE, NO, FI</div>