



Product Code	32511
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2529115
CND	C9005
EAN/UPC	8023279325119
PARAF	903176663
GTIN	08023279325119
GMDN	13755

Description

WAN PLUS double head pediatric - LATEX FREE - pink
with inner spring binaural, superior tubing quality and 2 pairs of soft colour coordinated ear tips in different sizes.

Features a chrome-plated zinc-alloy chestpiece, a screwon plastic retaining ring, outer spring binaural and colour matched moulded PVC tubing with non-chill ring.

Assure optimum acoustic performances.

Tubing length: 63 cm.

Pediatric stethoscope: diaphragm Ø 35 mm, bell Ø 30 mm.

Supplied in a cardboard box.

Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.