



Product Code	32513
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2529115
CND	C9005
EAN/UPC	8023279325133
PARAF	903176752
GTIN	08023279325133
GMDN	13755

Description

WAN double head pediatric - LATEX FREE - light blue

Features a chrome-plated zinc-alloy chestpiece, a screw-on plastic retaining ring, outer spring binaural and colour matched moulded PVC tubing with non-chill ring.
Assure optimum acoustic performances.

Chestpiece: Ø 35 mm, bell Ø 30 mm.

Tubing length: 63 cm.

Supplied in a cardboard box.

Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.