


# IMPERMEABLE, ABSORBENT, PRE-CUT MATS 90x90 cm



Product Code	25150
Unit of sale	box of 70 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2184860
CND	T030104
EAN/UPC	8023279251500

## Description

### ABSORBENT PRE-CUT MATS

- Size: 90 x 90 cm
- Absorbency: 1,700 ml/m<sup>2</sup>

### IMPERMEABLE ANTI-SLIP ABSORBENT BLUE COLOUR MAT FOR FLOORS

Pre-cut mats to be positioned on the floor of the operating room or near the sink in the pre-operating theatre for holding liquids that could fall on the floor.

Mats have high absorbent capacity and optimal speed absorbency.

The mats distribute the absorbed liquid by capillary method and maintain it at his internal structure avoiding the lifts on surface.

The special film at rear side blocks the mat on the floor and keeps it stopped without using adhesive matters, which could release glue traces on the floor.

They are perfectly adherent to the ground, avoiding risk of accidents.