

# FOLDED NAPKINS - 33x45 cm green



Product Code	25230
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2442606
CND	T0399
EAN/UPC	8023279252309
PARAF	987472608
GMDN	40511

## Description

FOLDED NAPKINS - pack of 500 pcs  
Napkins 2+1 (2 plies cellulose + 1 ply polyethylene)

- Napkins 33x45 cm - green

Pre-coated folded napkins, made of 2 plies of pure cellulose wadding laminated with 1 ply of polyethylene.  
High absorbency and waterproofness to assure high comfort.

Made in Italy.