


FOLDED NAPKINS - 33x45 cm light blue



Product Code	25231
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2442607
CND	T0399
EAN/UPC	8023279252316
PARAF	987472572
GMDN	40511

Description

FOLDED NAPKINS - pack of 500 pcs
Napkins 2+1 (2 plies cellulose + 1 ply polyethylene)

- Napkins 33x45 cm - light blue

Pre-coated folded napkins, made of 2 plies of pure cellulose wadding laminated with 1 ply of polyethylene.
High absorbency and waterproofness to assure high comfort.

Made in Italy.