



Product Code	25282
Unit of sale	box of 5 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2382098
CND	M04010102
PARAF	913254645
GTIN	04064035107703
GMDN	34864

Description

3M™ TEGADERM™ + PAD - STERILE - 9x10 cm - box of 5

3M™ code: 3586NP

Transparent sterile dressing with absorbent, non-stick pad. Breathable and waterproof polyurethane coating ensures total protection from external contamination and reduces dressing changes, allowing patients to bathe and shower. It guarantees high absorption capacity and optimal skin transpiration, accelerating wound healing.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, GR.