

PHARMADOCT ELASTIC COMPRESS. BANDAGES 2.5m x 6cm - carton of 6 boxes



Product Code	25337
Unit of sale	1 carton
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2253334
CND	M0303010203
EAN/UPC	8017990161088
PARAF	926497393

Description

PHARMADOCT: first aid, medication and wellness.

- Adhesive elastic compression bandage 2.5 m x 6 cm
- Pieces per box: 1
- Minimum order: 6 box

Graphic of PharmaDOCT line boxes are subject to changes.

All products are in blister with 7 languages package (GB, FR, IT, ES, DE, PT, GR).