

NYTRILE GLOVES EXTRA LONG - 400 mm - medium



Product Code	25405
Unit of sale	box of 50 pcs
Minimum order	1
Type	PPE
Class	III

UK-REP	Not required
EAN/UPC	8016629009050

Description

Nitrile gloves - extra long - medium - blue

Disposable nitrile examination gloves powder-free and latex free.
The micro textured surface assures a very good grip.
Suitable for hospital use.
The 400 mm length and resistance offer high level of protection.

Packed in dispenser of 50.
Non sterile.
Medical device Class I.
PPE cat III.
EN455-EN374-AQL<1.5