

GLOVE DISPENSER - triple - stainless steel



Product Code	25432
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279254327
PARAF	987469525

Description

Glove dispenser - triple - stainless steel

Hold 3 standard boxes of 100 gloves.
Made of AISI 304 stainless steel.
Autoclavable at 121 to 134 °C.
Can be easily hung to the wall and into ambulances.
Screw included.

Suitable for S - M - L box of gloves together
Size: 25.2x9.5xh 37.7 cm