

# "NOVA" LED TORCH - white



Product Code	25440
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
--------	----

EAN/UPC	8023279254402
---------	---------------

GMDN	12347
------	-------

## Description

Professional penlights, also suitable for diagnostic, with bulb or white LED with quick clip shifter mechanism or push button (25634). Argenta and Elegance made of anodized aluminium, Nova in topquality ABS plastic housing. Supplied with 2x1.5V AAA batteries.

Nova torch  
Light: LED  
Light intensity: 25,000 m  
Material: plastic  
Supplied with batteries: 2 x AAA  
Colour: white  
ON/OFF by: Pocket clip

Made in Germany.