

"ARGENTA" LED TORCH - metal - blue



| | |
|---------------|----------------|
| Product Code | 25442 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|----|
| UK-REP | No |
|--------|----|

| | |
|---------|---------------|
| EAN/UPC | 8023279254426 |
| GMDN | 12347 |

Description

Professional penlights, also suitable for diagnostic, with bulb or white LED with quick clip shifter mechanism or push button (25634). Argenta and Elegance made of anodized aluminium, Nova in topquality ABS plastic housing. Supplied with 2x1.5V AAA batteries.

Argenta torch
Light: LED
Light intensity: 25,000 m
Material: metal
Supplied with batteries: 2 x AAA
Colour: blue
ON/OFF by: Pocket clip

Made in Germany.