

SILICONE TUBE 3x5 mm - 1 mm thick



Product Code	25452
Unit of sale	roll of 100 m
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2207976
CND	V9099
EAN/UPC	8023279254525
PARAF	987473131
GTIN	08023279254525
GMDN	35917

Description

SILICONE TUBES for external medical use - silicone 100%

Resistant to hot air up to 200°C.
Infra-red vulcanization. Hours treatment at 200°C temperature in special oven. Anti-adherent, avoid incrustation and coagulation.
Medical Grade - Transparent - Odourless

Made in European Union according to: US PHARMACOPOEIA XXI EUROPEAN PHARMACOPEA DIN 58367 - DIN 58362 - DIN 13098- FDA 177.2600 - BGA XV part A and B - BGA IX part B

Made in Spain.

Technical Specifications

Raw material: HCR
Colour: translucent
Hardness: 60° ±5% Shore - DIN 53505
Density: 1.17 g/ccm - DIN 53479A
Stretch: 490/530% - DIN 53504S1
Tear resistance: 34-35 N/mm - ASTM D 624 B