



|               |                |
|---------------|----------------|
| Product Code  | 25505          |
| Unit of sale  | box of 50 pcs  |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | Yes           |
| RDM (NSIS) | 2510979       |
| CND        | V9001         |
| EAN/UPC    | 8023279255058 |
| PARAF      | 976329401     |
| GMDN       | 42461         |

## Description

### PLASTIC COLOURED TONGUE DEPRESSORS

Plastic tongue depressor animal shaped with bright colour.  
Each tongue depressor is individually packed in transparent pouch.  
Packed in box of 50 pieces.

Multilingual box: GB, FR, IT, ES, DE, PT, GR, Arabic