

WOOD TONGUE DEPRESSORS - not sterile



Product Code	25509
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2501482
CND	V9001
EAN/UPC	8023279255096
PARAF	950850608

Description

WOOD TONGUE DEPRESSOR - NON STERILE - box of 100

Wooden medical tongue depressors made of quality birch wood. Edges accurately finished.
Size: 150x18x1.6 mm

Multilingual box: GB, FR, IT, ES, PT, DE, GR, Arabic