

# WOOD TONGUE DEPRESSORS - not sterile (50 packs of 100)



Product Code	25510
Unit of sale	box of 5000 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2501482
CND	V9001
EAN/UPC	8023279255102
PARAF	903144261

## Description

WOOD TONGUE DEPRESSOR - NON STERILE - box of 5,000 (50 boxes of 100 pcs)

Wooden medical tongue depressors made of quality birch wood. Edges accurately finished.  
Size: 150x18x1.6 mm

Multilingual box: GB, FR, IT, ES, PT, DE, GR, Arabic