

PLASTIC TONGUE DEPRESSOR - adult/pediatric - non sterile



Product Code	25513
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	5605622433813
PARAF	987471087

Description

PLASTIC TONGUE DEPRESSORS ADULT/PEDIATRIC NON-STERILE

- Tongue depressor - box of 100

Translucent tongue depressors made of polystyrene with curved shape and smooth edges.
One end in adult size and one end in child size.
Latex and PVC free.

Curved shape allows perfect vision during examination avoiding that the hand of examiner obstruct the vision

- Child side 15 mm
- Adult side 20 mm
- Length 150 mm