



Product Code	25526
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2449888
CND	V9003
EAN/UPC	8023279255263
GTIN	08023279255263
GMDN	67045

Description

MILITARY TOURNIQUET - orange

Mainly used for stopping massive haemorrhage of limbs on field conditions or in normal pre-hospital emergency.

Disposable easy to use, one-handed application tourniquet improving the efficiency of hemostasis.

It can be self-applied or applied by a first responder, allowing a quick and easy regulation of the pressure needed.

Suitable for both upper and lower limbs.

Made of woven synthetic fabric with thermoplastic clasp and secure windlass.

Size: 95x4 cm.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.