



Product Code	25529
Unit of sale	box of 24 pcs
Minimum order	1
Type	Medical device
Class	I M

PARAF	989322247
GTIN	14018674454688

Description

CLASSIC THERMOMETERS display of 24 boxes

Classic therm - GB, FR, IT, ES, PT, DE, NL, NO, DK, FI, SE

Mercury-free glass thermometers developed by Geratherm® Medical AG, containing a patented liquid alloy consisting of gallium, indium and tin-called Galinstan. Galinstan is completely non-toxic and environment-friendly.

They are ideal for determining body temperature of people of all ages, suitable for rectal, oral or axillary measurement. Each thermometer is calibrated to guarantee measurement accuracy.

- measuring range: 35.5°C to 42°C
- sustainable compared to digital thermometers
- easy to disinfect thanks to the closed glass body, ideal for allergy sufferers
- reset tool included "Easy flip" shaking aid

Thermometers include a reset tool "easy flip" which makes it easier to lead back the measuring liquid of the thermometers by considerably less effort when shaking down.

Supplied with individual box protective case, reset tool and multilingual manual. Thermometers are sold in multilingual display of 24 pieces.

Made in Germany.