



| | |
|---------------|----------------|
| Product Code | 25532 |
| Unit of sale | box of 24 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I M |

| | |
|-------|----------------|
| PARAF | 989322211 |
| GTIN | 14018674454992 |

Description

CLASSIC+ THERMOMETERS display of 24 boxes

Classic+ therm - GB, FR, IT, ES, PT, DE, NL, NO, DK, FI, SE

Mercury-free glass thermometers developed by Geratherm® Medical AG, containing a patented liquid alloy consisting of gallium, indium and tin-called Galinstan. Galinstan is completely non-toxic and environment-friendly.

They are ideal for determining body temperature of people of all ages, suitable for rectal, oral or axillary measurement. Each thermometer is calibrated to guarantee measurement accuracy.

- measuring range: 35.5°C to 42°C
- sustainable compared to digital thermometers
- easy to disinfect thanks to the closed glass body, ideal for allergy sufferers
- reset tool included "Easy flip" shaking aid

Thermometers include a reset tool "easy flip" which makes it easier to lead back the measuring liquid of the thermometers by considerably less effort when shaking down.

Storage case with enlargement function

Special storage case with magnifying glass to help people read results more clearly

Supplied with individual box protective case, reset tool, storage case with magnifying lens and multilingual manual. Thermometers are sold in multilingual display of 24 pieces.

Made in Germany.