

# RI-THERMO SENSIOPRO NON CONTACT INFRARED THERMOMETER



Product Code	25540
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2877452
CND	V0301010202
EAN/UPC	4045396219953

## Description

### RI-THERMO SENSIOPRO NON-CONTACT INFRARED THERMOMETER

Contactless, accurate and fast thermometer, helping minimize cross contamination. Clinical grade forehead temperature measurements, calibrated to an oral body temperature reference. Optimized for professional.

#### Features:

- clinical grade reading in 1-2 seconds
- ISO standards 80601-2-56:2017
- safe and efficient
- infrared non-contact helps preventing cross contamination
- 3000 measurements on one set of batteries
- 30 measurement storage
- backlit LED screen
- single hand operation

Multilingual manual: GB, FR, IT, ES, DE, RU.

## Technical Specifications

Measurement range: Forehead: 22 - 44°C  
Accuracy: Forehead: 36 - 39°C ±0.2°C, < 36°C or > 39°C  
Operating Temperature: 10 to 40°C (50 to 104°F), 95% RH or less  
Resolution: 0.1°C / 0.1°F  
Measurement Unit: °C (Default) / °F  
Hi/Lo Indicator: 34 to 43°C (93,2 to 109°F)  
Power Source: Two 1.5 V AA Batteries (non-rechargeable)  
Battery Life: Approximately 3000 measurements  
Memory Capacity: Allows memory recall, 30 previous measurements  
Power saving: Automatic power off after 30 seconds  
Measuring distance: Approx. 3 - 7 cm  
Dimension / Weight: 150 (L) x 48.5 (W) x 55 (H) mm / 125.8 g (with battery)  
Standards: ASTM E1965-98 and EC directive 93 / 42 / EEC