



Product Code	25548
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1323879
CND	V03010199
EAN/UPC	8023279255485
UMDNS	14034
PARAF	973922495

Description

TEMP SITTER - baby wireless thermometer

Accurate and innovative thin thermometer.

Transmission to Temp Sitter App (iOS and Android).

To be placed under baby armpit.

Update frequency every 3 seconds.

High temp alarm can be set in the App.

Waterproof and material of safe quality.

Supplied with 10 adhesives to fix thermometer to baby armpit, manual (GB, FR, IT, ES).

Technical Specifications

Battery life: 3,000 hours, not replaceable

Measurement range: 25-45 °C (± 0.05 between 35-38.5 °C)

Size: 51x31.8x8.5 mm - Weight: 9.6 g

Switchable °C/°F