

GIMA BL1 WIDE SCREEN DIGITAL THERMOMETER °C - box



Product Code	25553
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1831018
CND	V03010102
EAN/UPC	8023279255539
UMDNS	14032
PARAF	973922584

Description

BL1 WIDE SCREEN DIGITAL THERMOMETER - box flexible tip, water-proof, backlight after measurement, wide di-splay (27 x 22 mm).

Battery 3 V button size. Low battery indicator.

Provide accurate and quick read-out. Safe and easy to use with clear and large display.

- display range 32.0°C- 42.9°C
- accuracy $\pm 0.1^\circ\text{C}$
- memory: last reading
- beeper function and auto shut-off
- plastic case
- multilanguage manual and box (GB, FR, IT, ES, PT, GR)