



Product Code	25556
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1202916
CND	V03010102
EAN/UPC	8023279255560
UMDNS	14032
PARAF	976329856
GMDN	14032

## Description

INSTANT DIGITAL THERMOMETER °C / °F switchable

- hang box, flexible tip, water-proof, fever alarm, large display

Provide accurate and quick read-out. Safe and easy to use with clear and large display.

- display range 32.0°C- 42.9°C
- accuracy  $\pm 0.1^{\circ}\text{C}$
- memory: last reading
- beeper function and auto shut-off
- plastic case
- multilanguage manual and box (GB, FR, IT, ES, PT, GR, PL, Arab)

Works with battery LR/4I - 1.5V