

# JUMBO 2 DIGITAL THERMOMETER °C - hang box



Product Code	25557
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1984763
CND	V03010102
EAN/UPC	8023279255577
UMDNS	14035
PARAF	981516305

## Description

FLEXI JUMBO 2 DIGITAL THERMOMETER °C - hang box, fl exible tip, water-proof resistant, fever alarm

Provide accurate and quick read-out. Safe and easy to use with clear and large display.

- display range 32.0°C- 42.9°C
- accuracy  $\pm 0.1^\circ\text{C}$
- memory: last reading
- beeper function and auto shut-off
- plastic case
- multilanguage manual and box (GB, FR, IT, ES, PT, GR, PL, Arab)

Works with battery LR/4I - 1.5V