

JUMBO DIGITAL THERMOMETER °C rectal/oral - hang box



Product Code	25558
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	930777
CND	V03010102
EAN/UPC	8023279255584
UMDNS	14-032
PARAF	927604494

Description

FLEXI JUMBO DIGITAL THERMOMETER °C / °F switchable

- hang box, flexible tip, water-proof, fever alarm, backlight jumbo display (27.5 x 12 mm size). Rectal/oral use.

Battery CR1632.

Provide accurate and quick read-out. Safe and easy to use with clear and large display.

- display range 32.0°C- 42.0°C
- accuracy $\pm 0.1^{\circ}\text{C}$
- memory: last reading
- beeper function and auto shut-off
- plastic case
- multilanguage manual and box (GB, FR, IT, ES, PT, GR)