



Product Code	25559
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1967553
CND	V03010102
EAN/UPC	802327925591
UMDNS	14-032
PARAF	980908204

Description

GIMA DIGITAL THERMOMETER °C - box - 30 seconds
Sold in a counter display of 36 pcs. (GB - IT)

Provide accurate and quick read-out. Safe and easy to use with clear and large display.

- display range 32.0°C- 42.9°C
- accuracy $\pm 0.1^{\circ}\text{C}$
- memory: last reading
- beeper function and auto shut-off
- plastic case
- multilanguage manual and box (GB, FR, IT, ES, PT, GR)

Works with battery LR/4I - 1.5V