

# GIMA DIGITAL THERMOMETER °C - hang box - GB,FR,IT,ES,PT,DE,GR,Arabic



Product Code	25560
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	17645
CND	V03010102
EAN/UPC	8023279255607
PARAF	905898363
GTIN	08023279255607
GMDN	14035

## Description

GIMA DIGITAL THERMOMETER °C - hang box - 60 seconds

Provide accurate and quick read-out. Safe and easy to use with clear and large display.

- display range 32.0°C- 42.9°C
- accuracy  $\pm 0.1^{\circ}\text{C}$
- memory: last reading
- beeper function and auto shut-off
- plastic case
- multilanguage manual and box (GB, FR, IT, ES, PT, DE, DK, LT, GR, DE, Arabic)

Works with battery LR/4I - 1.5V