



Product Code	25561
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	17645
CND	V03010102
EAN/UPC	8023279255614
PARAF	987472990
GTIN	08023279255614
GMDN	14035

Description

GIMA DIGITAL THERMOMETER °F - hang box - 60 seconds (89,6-107,6 °F)

Provide accurate and quick read-out. Safe and easy to use with clear and large display.

- accuracy $\pm 2,0^{\circ}\text{F}$
- memory: last reading
- beeper function and auto shut-off
- plastic case
- multilanguage manual and box (GB, FR, IT, ES, PT, DE, GR, Arabic)

Works with battery LR/4I - 1.5V