



Product Code	25575
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	894174
CND	V03010102
EAN/UPC	8023279255751
PARAF	977810720

## Description

### PROFESSIONAL INFRARED EAR THERMOMETER

- special design for hospital use
- durable probe for cover replacement up to 100,000 times
- probe cover ejector and probe cover detector
- one-second response time: fast and easy
- neck strap and 1 probe cover included
- fever alarm
- selectable °C/°F

Needs probe covers for use

Large LCD display

Multilingual box and manual: GB, FR, IT, ES.

## Technical Specifications

Measurement range: 34~42.2 °C (93.2-108 °F)

Accuracy: ±0.2 °C (0.4 °F) during 35.5~42 °C (95.9~107.6 °F), ±0.3 °C (0.5 °F) for other range

Memory: 25 sets

Buzzer: Yes

Battery: CR2032 1 (Gima code 32212)

Dimension - weight: 144.82x39.84x48.50 mm - 68.8 g