



Product Code	25579
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	912714
CND	V03010102
EAN/UPC	8023279255799
UMDNS	17-887
PARAF	924994686

Description

POCKET EAR THERMOMETER

Small, easy to clean, cover free design ear thermometer. Home use.

Easy to read 3 digit LCD display 27.9x13 mm

Multilanguage box and manual: GB, FR, IT, ES.

Technical Specifications

Measurement Range: 34-44 °C (93.2-111.2 °F) - ± 0.2 °C (± 0.4 °F)

Response time: 1-2 seconds

Battery: Lithium 3 V (CR2032) - life approx 4,000 measurements

Size - weight: 114x35x40 mm - 35 g (including batteries)