



Product Code	25580
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes

RDM (NSIS)	237591
CND	V03010102
EAN/UPC	8023279255805
PARAF	903176358
GTIN	08023279255805
GMDN	17887

Description

IR EAR THERMOMETER

- fully meet "ASTM" standard and CE EN 12470-5 (only few brands passed it) in measurement accuracy (35.5 °C to 42 °C ± 0.2 °C, others ± 0.3 °C).
- unique and handy design
- one-second measurement, fast and easy
- 9 memory for temperature comparison and evaluation, fever alarm, easy reading LCD display
- self-diagnostic and infrared measurement

Multilingual box and manual: GB, FR, IT, ES, PT, DE, GR, Arabic* * Only box

Technical Specifications

Measurement Range: 34.0 °C~42.2 °C (93.2 °F~108.0 °F)
 Response time: 1 second
 Time between measurements: > 5 sec
 Display Model °C/°F
 Battery Type: CR2032 (32212)
 Dimension (LxWxH) - Weight: 14x3.8x3 cm - 70 g
 Norms: IEC60601-1, IEC60601-1-2 (EMC), ASTM E1965-98