

AEON A200 NON CONTACT INFRARED THERMOMETER GB,FR,IT,ES



Product Code	25582
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1937072
CND	V03010102
EAN/UPC	8023279255829
UMDNS	14036
PARAF	980515353
GMDN	17888

Description

A200 INFRARED NON CONTACT THERMOMETER - GB, FR, IT, ES

Multifunctional infrared non contact thermometer for body and objects.

Effective distance 2-3 cm.

- result in 1 second
- LCD backlight display
- quality sensor for high accuracy measure
- °C and °F are easily convertible
- 3 colours screen red, orange and green indicate temperature level
- automatically power off after 30s if left it on still
- 25 groups memories
- audible alarm
- low power warning

Supplied with multilingual box and manual (GB, FR, IT, ES or IT, GR, PL, RO) and batteries.

Body or objects switch.

- Green backlight < 37.5°C (99.5 °F)
- Orange backlight = 37.5°C = 37.9°C (=99.5°F =100.2°F)
- Red backlight = 38.0°C (=100.4 °F)

Technical Specifications

Measurement Range:

- Body 32 - 42.9°C (89.6 - 109.3 °F) ±0.2 °C
- Objects 0-99°C (32.0-211.82 °F)

Response time: 1 second

Memory: 25 groups

Display Mode: °C/°F

Battery Type: 2 AAA 1.5V

Dimension (LxWxH) - Weight: 168x54x44 mm - 75 g