

GIMATEMP NO CONTACT INFRARED THERMOMETER GB,FR,IT,ES



Product Code	25583
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1951667
CND	V03010102
EAN/UPC	8023279255836
PARAF	980505604

Description

GIMATEMP NON CONTACT THERMOMETER GB, FR, IT, ES

Multifunctional infrared non contact thermometer with 4 applications: body, liquid, food, ambient temperature.

- high accuracy and quick measurement in =2-3 sec
- 3 colour backlight alarm (red, yellow, green) for checking the temperature status
- measuring distance 1-5 cm
- room temperature is always displayed.
- 50 memories

Technical Specifications

Measurement method: Temporal artery scan or forehead spot

Display model - resolution: °C/°F - 0.1°C/0.1°F

Operating mode: Adjusted mode (Body mode)

Direct mode (Surface mode)

Output range:

- Body mode: 34.0~43.0°C/93.2-109.4°F

- Surface mode: 0~100°C/32-212°F

Accuracy: Body mode: 35.0~42.0°C ±0.2°C / 95.0-107.6°F ±0.4°F

±0.3°C/±0.5°F for other range

Surface mode: ±2°C/±3.6°F

Auto power off time: =18 s

Measuring time: 2-3 s

Measuring distance: 1-5 cm (0.4-2 inch)

Batteries (included): 1.5 V (AAA) alkaline batteryX2 (IEC type LR03)

Size - Weight: 122x96x42 mm - 80 g