



Product Code	25588
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279255881

## Description

### VET RECTAL THERMOMETER

Designed to be ergonomically correct for animals.  
Setting for pets, pigs, cattle and other animals.

#### Features:

- instant temperature in 6 seconds
- beeps indicator
- LCD display and memory
- automatic shut-off after 15 seconds
- display temperature in both Fahrenheit or Celsius scale
- reliable, safe and easy to use

Multilingual box and manual: GB, FR, IT, ES, PT, GR.

## Technical Specifications

Measurement range: 32~44°C (89.6-111.2 °F)  
Accuracy: ± 0.1 °C (35~39°C)  
Battery life is approx. 5,000 readouts  
Dimension: 190x110x35 mm  
Battery: 2x1.5V AAA  
Measurement time: 6 sec.