



| | |
|---------------|----------------|
| Product Code | 25600 |
| Unit of sale | bag of 10 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I M |

| | |
|------------|---------------|
| RDM (NSIS) | 699227 |
| CND | V03010199 |
| EAN/UPC | 8008687000223 |
| PARAF | 920010941 |

Description

FEVERCONTROL - bag of 10 blister

Reusable, non toxic, mercury free, safe, unbreakable, reliable, long lasting, easy to read, liquid crystal thermometers. They indicate temperature in 15 sec with graduation 0.5°C. Medical Devices.

Made in Italy.

Blisters in GB, IT