



Product Code	25601
Unit of sale	bag of 10 pcs
Minimum order	1
Type	Medical device
Class	I M

RDM (NSIS)	1243453
CND	V03010199
EAN/UPC	8008687000117
PARAF	906985561

Description

48h FEVER MONITOR - bag of 10 blister of 3 pcs.

Innovative "continual" fever indicator for 48h monitoring, for babies.
Disposable, non toxic, mercury free, safe, unbreakable, reliable, easy to read, liquid crystal thermometers.
It indicates temperature in 30 seconds.
Medical Device.

Made in Italy.

Blister in GB, IT