

OMRON FLEX TEMP SMART DIGITAL THERMOMETER - MC 343F-E



Product Code	25604
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2331939
CND	V0301010201
EAN/UPC	4015672113602

Description

FLEX TEMP SMART THERMOMETER - MC 343F-E

Technical Specifications

Measurement Type Oral, Armpit, Rectal
Flexible tip
Beeper (signals when a measurement is complete)
Memory last reading
°C/°F selection
Large, 3 digit display
Water resistant
Automatic shut off