

DIGITAL BASAL THERMOMETER °C - plastic case - 2 decimals



Product Code	25608
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	237594
CND	V03010102
EAN/UPC	8023279256086
PARAF	913758619
GTIN	08023279256086
GMDN	14035

Description

Digital basal thermometer.

With 4 digits, 2 decimals, for fertility control.
Water-proof.

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR, Arab.