

OMRON GENTLE TEMP 520 EAR THERMOMETER - MC-520-E



Product Code	25610
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A
UK-REP	No
CH-REP	No
RDM (NSIS)	1660080
CND	V03010102
EAN/UPC	4015672107793

Description

OMRON GENTLE TEMP® 520 EAR THERMOMETER

Gentle Temp® 520 infrared thermometer offers comfortable, safe, accurate and quick temperature measurement from the eardrum. It offers advanced features: one-second measurement and a memory for up to nine readings to track temperature changes over a period of time and monitor trends.

It is intended for home use.

Large LCD display.

Switchable reading scale °C/°F.

It can measure consecutively up to 3 times.

Supplied with 1 battery, 21 probe covers, probe cap, connection ring and multilingual instructions (GB, FR, IT, ES, DE, NL, RU, TR, Arabic).

Technical Specifications

Measurement range: 34.0-42.2°C (93.2-108.0°F)

Accuracy: ± 0.2°C (± 0.4°F)

Measurement time: 1 second

Temperature mode: °C/°F switchable

Auto power off time: 1 min

Memory: 9 readings

Display: LCD display

Battery: 3.0V DC, 1 CR2032 (Gima code 32212) lithium button battery

Battery life: approx. 2,500 measurements

Size - weight: 161x36x56 mm - 85g

Classification - protection: Type BF - IP22