

OMRON GENTLE TEMP 720 NON CONTACT THERMOMETER - MC-720-E



Product Code	25611
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A
UK-REP	No
CH-REP	No
RDM (NSIS)	1462417
CND	V03010102
EAN/UPC	4015672108561

Description

OMRON GENTLE TEMP® 720 NON CONTACT THERMOMETER

Gentle Temp® 720 is a superfast accurate infrared non contact thermometer. It is intended to measure body, surface and room temperature. It can memorise up to 25 readings, track temperature changes over a period of time and monitor trends. In addition to the forehead mode, it can also be switched to surface measurement mode, which shows the actual and unadjusted surface temperature.

It helps to determine if the object temperature (for example a milk bottle or food dispenser) is suitable for the patient or the baby. Room temperature is always displayed. It is intended for home use. Reading scale °C/°F is switchable. Supplied with 1 battery, probe cap and multilingual instructions manual: GB, FR, IT, ES, DE, NL, RU, TR, Arabic.

3 APPLICATIONS: FOREHEAD, SURFACE, ROOM

Technical Specifications

Measurement type: contactless 3 in 1 (forehead, surface, room)
Measurement range/accuracy:
- forehead 34.0 - 42.2°C (93.2-108.0°F) / ± 0.2°C
- surface -22.0 - 80.0°C (-7.6° - 176.0°F) / ± 0.3°C
- room 10 - 40°C (50°-104°F) / ± 2°C
Auto power off time: 1 min
Measuring time: 1 second
Measuring distance: 1-3 cm (0.4-1.2 inch)
Temperature mode: °C/°F
Memory: 25 readings
Display: backlight
Battery: 3.0V DC, 1 CR2032 (code 32212) lithium button battery
battery life: approx. 2,500 measurements
Classification - protection: Type BF - IP22
Size - weight: 155x39x45 mm - 90g