



Product Code	25615
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2306429
CND	V0301010202
EAN/UPC	8023279256154
GMDN	63711

Description

BLUETOOTH INFRARED FOREHEAD AND EAR THERMOMETER

Multifunctional infrared thermometer used to read body (ear and forehead), room and object temperature.

- clinically accurate, reliable and quick
- professional and home use
- patent technology to separately measure adults and children
- four colour LED backlight
- 10 memory records
- switchable between °C and °F
- audio and visual fever warning
- magnetic probe cover to automatically switch between forehead and ear mode

Bluetooth connection to mobile devices, downloading free Vihealth App.

Supplied with pouch, 2 AAA batteries, multilingual manual: GB, FR, IT, ES, DE.

Technical Specifications

Measurement type: 3 in 1 (forehead, ear, object, liquid)
 Measurement range/accuracy: forehead & ear 32.0-42.9°C ± 0.2°C / object 0 -100°C / ± 1.0°C
 Auto power off time: 1 min
 Measuring time: 1 second
 Temperature mode: °C/°F switchable
 Memory: 10 groups
 Display: LCD, four-colour LED display (white, green, orange, red)
 Resolution: 0.1 °C/°F
 Battery: 3V DC 2xAAA batteries, can be used more than 3000 times
 Size - weight: 143x35x41 mm - 80g (with batteries)
 Classification - protection: Type BF - IP22