



Product Code	25641
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2671671
CND	L1799
EAN/UPC	8023279256413
GTIN	08023279256413
GMDN	38832

Description

Spirit torch
 Light: LED
 Material: chrome plated brass
 Supplied with batteries: 2 x AAA
 Colour: silver
 Packaging: box
 Multilingual instructions: GB, FR , IT , ES
 ON/OFF by: Pocket clip

With pupil gauge and dual slide ruler.